


## Reduced Representation Bisulfite-Seq – A Brief Guide to RRBS

### What is RRBS?

Typically, RRBS samples are generated by digesting genomic DNA with the restriction endonuclease *MspI*. This is followed by end-repair, A-tailing, adapter ligation and finally bisulfite conversion. Often, the library is also size-selected for fragments between 40 and 220bp in length. This fragment size has been shown to be plentiful in the sample and yield information on the vast majority of CpG islands (CGIs) in the human or mouse genome. Fig. 1 shows that quite a few *MspI*-*MspI* fragments (generated in silico for the mouse genome) are even shorter than 40bp. Since the size selection process is not as good as it is in theory, often a sizeable number of fragments below 40bp can end up in the RRBS library.




The fairly small fragment size of RRBS fragments can become a potential problem especially for sequencing reads with high read length (e.g. > 75bp or >100bp). If the read length is



the same number of reads as both paired-end reads added together are more likely to yield more genuine methylation information, as long as the read length is long enough to allow for a fair mapping efficiency for single-end reads (40-50bp reads are probably long enough to get mapping efficiencies in the range of 60 to 70%).

## Other read length effects

Current sequencing on the Illumina platform often produces data whose quality deteriorates towards later cycles. We have received feedback from numerous sources, or downloaded data from public archives, which looks similar to this read quality profile:



Up to a read length of around 60-70bp, the basecall quality of these reads is excellent ( $>$  Phred 30). After that, however, Phred scores tend to drop dramatically in a fairly large number of sequences, which means the rates at which bases are called erroneously increases. Base call errors in reads can result in reads not being aligned at all (reduced mapping efficiency), incorrect methylation calls or, in the worst case, mis-alignments (which will most likely also generate incorrect methylation calls). Some of these aspects have been recently discussed in the following review: Krueger F., Kreck B. et al., DNA Methylation Analysis Using Short Bisulfite Read Data, 2012).


### In a nutshell:

RRBS reads suffer disproportionately from problems associated with long read lengths because they are, compared to other -Seq applications, size-selected for rather short fragment sizes. In the following sections, I will discuss some further aspects that need to be considered when analysing RRBS samples, and I will introduce our way of dealing with all sorts of read length related problems or experimentally introduced biases for RRBS libraries.

## Potential biases for directional RRBS libraries (OT, OB strands only)

Directional sequencing, the most wide-spread way of (RR)BS, only ever sequences reads originating from the original top (OT) or original bottom (OB) strands. For simplicity, I have refrained from drawing out all four bisulfite DNA strands for the illustration below.

The sequential steps of RRBS are:



**Cytosines in blue** retain the original genomic methylation state, whereas **cytosines in red** are introduced experimentally during the fragment end-repair reaction. This can be accomplished with either unmethylated or methylated cytosines, the trend seems to be that unmethylated cytosines are being used primarily now.

After the adapters are attached, the sequences are treated with sodium bisulfite, which converts unmethylated cytosines into thymines. Thus, the first three bases of (almost) all RRBS reads are either **CGG** or **TGG**, depending on their genomic methylation state. This applies to reads from both the OT and OB strand, and as nearly all reads in a directional RRBS experiment start with one of these two options, every read provides information on at least one CpG right in the start.






After bisulfite conversion, the first three bases of non-directional RRBS reads that originated from the OT or OB strands will also be either **CGG** or **TGG**, depending on their genomic methylation state. In addition, however, non-directional libraries may contain reads which originated from the CTOT or CTOB strands. These reads will have **CAA** or **CGA** at the start, depending on whether unmethylated or methylated cytosines were used for the end-repair reaction, respectively. (Theoretically, the sequence could also be **CAG** or **CGG**, but this would assume that a C in CHH context was methylated on the other strand, and this is arguably very rarely the case for CpG-rich sequences. For simplicity it is therefore left out here). In either case, the second base would incur a methylation call of a base that does no longer reflect the genomic methylation state, which is illustrated below.


For non-directional libraries, one can discriminate the following four cases:

**A) The read length is shorter than the MspI fragment, OT or OB alignments**



In this case, the entire read can be used for alignments and methylation calls. The first position resembles the true genomic methylation state (which can be C or T).

**B) The read length is shorter than the MspI fragment, CTOT or CTOB alignments**



In this case, the read will start with **CAA** or **CGA** (the filled-in and the other cytosine in CHG context on the opposing strand are expected to be fully bisulfite converted), whereby the position marked in **RED** infers the methylation state from a cytosine that was experimentally introduced. The positions in **BLUE** carries genomic methylation information. As a consequence, methylation information of the second base would bias the results depending on the methylation state of the cytosine used for end-repair and needs to be excluded from methylation analysis.






To make this calculation more reliable one can vary the required overlaps with the adapter sequence, e.g. requiring 3, 5, 7bp etc., before the non-conversion rate of the filled-in position is determined.

## Adaptive quality and adapter trimming with Trim Galore!

We have tried to implement a method to rid RRBS libraries (or potentially all kinds of sequencing datasets) of potential problems in one convenient process. For this we have developed a wrapper script (`trim_galore`) that makes use of the publically available adapter trimming tool `Cutadapt` and `FastQC` for optional quality control once the trimming process has completed.

Even though `trim_galore` works for any (base space) high throughput dataset (e.g. downloaded from the SRA) this section describes its use mainly with respect to RRBS libraries.


- In the first step, low-quality base calls are trimmed off from the 3' end of the reads before adapter removal. This efficiently removes poor quality portions of the reads. Here is an example of a dataset downloaded from the SRA which was trimmed with a Phred score threshold of 20 (data set DRR001650\_1 from Kobayashi et al., 2012).



before trimming

after trimming

In the next step, `Cutadapt` finds and removes adapter sequences from the 3' end of reads. If no sequence was supplied it will use the first 13 bp of the standard Illumina paired-end adapters ('AGATCGGAAGAGC'), which recognises and removes adapters from most standard libraries. To control the stringency of the adapter removal process one gets to specify the minimum number of required overlap with the adapter sequence, else it will default to 1. This default setting is extremely stringent, i.e. an overlap with the adapter sequence of even a single bp is spotted and removed. This may appear unnecessarily harsh, however, as a reminder adapter contamination may in a bisulfite-Seq setting lead to mis-alignments and hence incorrect methylation calls, or result in the removal of the sequence as a whole because of too many mismatches in the alignment process. Tolerating adapter contamination is most likely detrimental to the results, but we realize that this process may in some cases also remove some genuine genomic sequence. It is unlikely that the removed bits of sequence would have been involved in methylation calling anyway (since only the 4<sup>th</sup> and 5<sup>th</sup> adapter base would possibly be involved in methylation calls (for directional libraries that is)), however, it is quite likely that true adapter contamination – irrespective of its length – would be detrimental for the alignment or methylation call process, or both.



This example (same dataset as above) shows the dramatic effect of adapter contamination on the base composition of the analysed library, e.g. the C content rises from ~1% at the start of reads to around 22% (!) towards the end of reads. Adapter trimming with `Cutadapt` gets rid of most signs of adapter contamination efficiently. Note that the sharp decrease of A at the last position is a result of removing the adapter sequence very stringently, i.e. even a single trailing A at the end is removed.

- `trim_galore` also has an `'--rrbs'` option, which identifies sequences that were adapter-trimmed and removes another 2 bp from their 3' end. This is to avoid that the filled-in cytosine position (marked in **RED** in the above examples) close to the second `MspI` site in a sequence is used for methylation calls. Sequences which were merely trimmed because of poor quality will not be shortened any further.
- `trim_galore` also has a `'--non_directional'` option, which will screen adapter-trimmed sequences for the presence of either `CAA` or `CGA` at the start of sequences and

clip off the first 2 bases if found. If CAA or CGA are found at the start, no bases will be trimmed off from the 3' end even if the sequence had some contaminating adapter sequence removed (in this case the sequence read likely originated from either the CTOT or CTOB strand).

- Lastly, since quality and/or adapter trimming may result in very short sequences (sometimes as short as 0 bp), `trim_galore` can filter trimmed reads based on their sequence length (default: 20 bp). This is to reduce the size of the output file and to avoid crashes of alignment programs which require sequences with a certain minimum length.

**Paired-end reads:** Most alignment programs require both paired-end `FastQ` files to contain sequence pairs in exactly the same order. To not bring this sequence-by-sequence order out of sync, `Trim Galore` trims paired-end `FastQ` files sequentially but does not discard too short sequences straight away. Instead, once both files have been trimmed it performs a validation of both paired-end files at the same time where it can discard entire sequence pairs if one of the two sequences has become too short during the trimming. This process retains the sequence-by-sequence order of paired-end files, and results in the final output files of paired-end trimming, termed `_val_1.fq` and `_val_2.fq` (for validated). `Trim Galore` also allows you to keep singleton reads if they fulfil certain length criteria (see the `Trim Galore` help for more details).

Applying these steps to both self-generated and downloaded data can ensure that you really only use the high quality portion of the data for alignments and further downstream analyses and conclusions.